

 WINGDERM

# Presentation For RENUVA

Aesthetics & Technology, Easy to Achieve



# CONTENTS



1 RENUVA™

2 Cyroshot™

3 Global Marketing

4 Operation



**01**


**RENUVA™**

# RENUVA™ Introduce




## RENUVA™

- 1550nm Erbuim fiber laser
- Medical CE & FDA cleared
- The laser supply with long service life and high reliability
- Use for Skin Rejuvenation, Scar, Wrinkle, Hair Regrowth...




# RENUVA™ Introduce

RENUVA™ technical video




# RENUVA™ Introduce

- High Speed Pattern Generator
- Super small spot size, high speed, high precision scanning
- Scanning speed up to 100mm/s, precision up to  $\pm 0.1\text{mm}$
- Add random scan function, to reduce the pain, Faster recovery




# RENUVA™ Introduce



- Five Tips, Add newest roller tips;
- Automatic identification different tips
- Convenient for treatment in different areas
- Specially designed scalp treatment head, treatment scalp more fit
- Automatically match the treatment parameters according to the treatment head
- **Compare with Fraxel 1550nm**



# RENUVA™ Introduce



- Built-in Storage Box For Tips
- Convenient to store
- Rubber tips slot, prevent tips surface from wearing



## WingCloud

- IOT(Internet Of Things) Makes More Possibilities
- Economical Service Solutions
- Flexible Marketing
- More Sales opportunities


<http://cloud.wingderm.cn/Admin/Login/index>

IOT(Internet Of Things) Function;

1. New sales function developed for distributors, rent and control the device with Wingcloud.
2. Each distributor has an ID which is used to control and manage all the listed devices..
3. Wingcloud information, as below:

Device SN number; devices position; devices use time; device shot point; recharge management; device situation; and device failure checking(failure checking will launch end of 2021)

Order	SN	Name	Administrator	Location	Version	Configuration	Working Time (Hours)	Total usage	Balance	Status	Remarks	Remote operation	Operation		
1	L0200910004	808 laser hair removal instrument	Tsvetan Ivanov	Bulgaria	LD.A-1.1.1...	WaveLength	0.64hours	791points	9279points	Shutdown		Lock	Remote Control	Reset	History
2	L0200810002	808 laser hair removal instrument	Ali Ghazafar	South Africa	LD.A-1.1.1...	1000W, 3 Wa...	Hours	0points	1000points	Shutdown		Lock	Remote Control	Reset	History
3	L0200810001	808 laser hair removal instrument	Ali Ghazafar	South Africa	LD.A-1.1.1...	1000W, 3 Wa...	Hours	0points	1000points	Shutdown		Lock	Remote Control	Reset	History
4	L0200710019	808 laser hair removal instrument	xavier	France	LD.A-1.1.1...	1000W, Th...	75.02hours	23649points	47633points	Shutdown	IREM COFF...	Lock	Remote Control	Reset	History
5	L0200510007	808 laser hair removal instrument	xavier	France	N/A	1000W/3 Wa...	Hours	0points	0points	Shutdown	1000W/31 ...	Lock	Remote Control	Reset	History
6	L0200210008	808 laser hair removal instrument	xavier	France											
7	L0190910016	808 laser hair removal instrument	xavier	France											
8	L0181240001	808 laser hair removal instrument	test_watgion	China											
9	L0200510001	808 laser hair removal instrument	ERD	China											
10	L0191010001	808 laser hair removal instrument	test_xzy	China											





02

CryoShot™

# Cyroshot Introduce



## CRYOSHOT™



- **Ultra-low temperature**
- High air volume
- Flexible adjustment bracket
- **Connection with RENUVA,**  
Could be also match with IPL, CO2, or other laser

# Cryoshot™ Introduce

## Smart Design





Exquisite and precise operation




# Cryoshot™ Introduce

## Adsorbable Nozzle Match with RENUVA




# Cryoshot™ Introduce



- **Ultra-low temperature**
- Core temperature is as low as -40 °C
- Cooling at -20°C
- **Use Hitachi's high-efficiency compressor**
- Honey comb heat exchange system
- Better skin protection and analgesic effect

# Cryoshot™ Introduce



- **Automatic defrost**
- **High air volume**
- Up to 800L/min high air volume
- Ensure 8 hours continuous work
- **Compare with Cyro 6**





**03**

**Global Marketing**




# Global Marketing

## RENUVA™



© Coline Marslen



# Global Marketing

## RENUVA™ Feedbacks



# Global Marketing

## RENUVA™ Feedbacks



# Global Marketing


UK distributor



# Global Marketing

RENUVA™

First online Show 2021






**04**

**Operation**


# Operation

## Scan Mode

Scrambled Mode: Laser shooting is uncontrolled with no repetition



Random Mode: Laser shooting is controlled with repeatable



# Operation

## Built-in Parameters



### Parameter

- Energy: 1mJ ~ 70mJ
- Shots Density: 100 Shots/cm<sup>2</sup> ~ 500 Shots/cm<sup>2</sup>



# Thanks



**WINGDERM**

[www.wingderm.com](http://www.wingderm.com)